




Clinical Data



C-STEM®
TOTAL HIP SYSTEM



98.96%
survivorship at 14 years¹

- The original C-STEM® Implant established a tapered slip stem that not only worked with the bone cement mantle, but was also designed to produce positive proximal bone remodelling due to more anatomic stress distribution.²
- C-STEM & ELITE PLUS OGEE® Cemented Cup are shown to have a 97,37% survival rate at 10 years in over 5,100 cases³
- In their patient series follow up, Sochart et al. (2009) noted "there was no evidence of negative bone modelling", which is an issue seen with other cemented stems.⁴

PINNACLE®
HIP SOLUTIONS



99.2%
survivorship at 10 years⁵

- With nearly 15 years clinical use the PINNACLE® Acetabular Cup System is DePuy Synthes' premium Acetabular System.
- PINNACLE Acetabular Cup System offers surgeons the ability to address existing pathologies whilst choosing from MARATHON® or ALTRX® cross linked polyethylene or CERAMAX™ Ceramic Liners.
- 10 years data in combination with C-STEM AMT Stem in the NJR with a survival rate of 96.21%³



Launched in 2005* the C-STEM® AMT has been provided to patients in over 100,000 total hip replacements Worldwide with 50,000 in the UK alone.⁶

* C-STEM AMT standard and high offset stems launched in 2005.



In 2017 the C-STEM AMT Total Hip System was awarded an ODEP 10A* rating by the Orthopaedic Data Evaluation Panel.⁷

98.64%³

C-STEM AMT and MARATHON Cemented Cup is shown to have a 98.64% survival rate at 7 years in over 8,100 cases.³ In combination with the PINNACLE Cementless Cup it is shown to have a 96.21% survival rate at 10 years in over 10,000 cases.³

Combination	Implantations	Year 1	Year 3	Year 5	Year 7	Year 10
C-STEM AMT Cemented Stem/ MARATHON	8,111	0.43 (0.30-0.61)	0.94 (0.71-1.24)	1.23 (0.89-1.69)	<i>1.36 (0.96-1.92)</i>	-
C-STEM AMT Cemented Stem/ PINNACLE	10,027	0.70 (0.55-0.89)	1.22 (0.99-1.49)	1.75 (1.42-2.16)	2.23 (1.76-2.83)	3.79 (2.81-5.12)

National Joint Registry for England, Wales, Northern Ireland and the Isle of Man. Revision rates (all-cause) for main hip stem and cup combinations (95% confidence intervals).

Note: Blue italics signify that fewer than 250 cases remained at risk at these time points

94.3%⁸

C-STEM AMT, used in combination with the PINNACLE cementless cup is shown to have a 94.3% survival rate at 10 years in more than 3,000 cases.⁸

Combination	Implantations	Year 1	Year 3	Year 5	Year 10
C-STEM AMT Cemented Stem/ PINNACLE	3,025	1.3 (0.9, 1.8)	2.4 (1.8, 3.1)	3.1 (2.4, 4.1)	5.7 (3.6, 8.8)

Extracted from Table HT13. Australian Orthopaedic Association National Joint Replacement Registry. Cumulative Percent Revision of Primary Total Conventional Hip Replacement with Hybrid Fixation (95% confidence intervals).



96.9%
survivorship at 99 months⁹

- The C-STEM AMT builds upon the clinical heritage of the original C-STEM Hip and its specific triple tapered, highly polished design. The functional intra-medullary geometry of original C-STEM Stem has been maintained whilst the extra-medullary geometry has been enhanced with a raised lateral shoulder, designed to improve visualization during insertion, and the addition of the 12/14 ARTICUL/EZE™ mini-taper.
- Berstock et al (2014) followed 415 C-STEM AMT Hips in 386 patients for 5-8 years. There were only 2 femoral stem revisions, both for infection, and the Kaplan-Meier survivorship at 99 months was 96.9%.⁹
- Flatoy et al (2015) followed 70 hips using C-STEM AMT and MARATHON Cemented Cups utilising RSA analysis. At 2 years follow up, there were no revisions and authors concluded that the C-Stem AMT Stem displayed a pattern of migration and bone remodeling that predicts good clinical performance.¹⁰

MARATHON®
CROSS-LINKED POLYETHYLENE



No failures
at 5 years¹¹

- Introduced in 2008, MARATHON Cemented Cups utilise the cross linked polyethylene successfully used for over 10 years as a modular liner in the DURALOC® and PINNACLE Acetabular Systems.^{5,12}
- MARATHON® Cross-Linked Polyethylene has been shown to be significantly more resistant to wear compared to non-crosslinked polyethylene - in 41 cases, Stone et al (2014) achieved 100% survivorship at 5.2 years with 83% reduction in wear versus ENDURON™ POLYETHYLENE.¹¹
- The fixation features and bearing surface geometry of the MARATHON Cemented Cup are the same as the LPW version of the CHARNLEY® OGEE, CHARNLEY & ELITE PLUS OGEE CUPS.

References

1. Wroblewski BM, et al. The C-Stem in Clinical Practice, Fifteen-Year Follow-Up of a Triple Tapered Polished Cemented Stem. *J.Arthroplasty*. 2013; 28: 1367-1371
2. Wroblewski BM, et al. Triple Tapered Polished Cemented Stem in THA. *J Arthrop*, 2001.
3. National Joint Registry for England, Wales, Northern Ireland and the Isle of Man, 15th Annual Report, 2018. Table 3.9. Available from www.njrreports.org.uk.
4. Sochart D,H. et al. Results of the C-STEM polished triple tapered femoral component. *J.BJS Br* 2010; 92-B:SUPP I 99
5. Callaghan J,J. et al. Fixation and Wear With a Contemporary Acetabular Component and Cross-Linked Polyethylene at Minimum 10-Year Follow-Up. *JOA* 2014; 29: 1961-69
6. C-Stem Sales Year 2005 to Year 2018 data on file DePuy Synthes.
7. Orthopaedic Data Evaluation Panel. ODEP product ratings. Latest ODEP ratings can be found at www.odep.org.uk [Accessed 01/11/2018].
8. Australian Orthopaedic Association National Joint Replacement Registry. Annual Report. Adelaide: AOA; 2018. Available from URL: <https://aoanjrr.sahmri.com/annual-reports-2018>. Extracted from Table HT13.

Table HT13 Cumulative Percent Revision of Primary Total Conventional Hip Replacement with Hybrid Fixation

Femoral Component	Acetabular Component	N Revised	N Total	1 Yr	3 Yrs	5 Yrs	7 Yrs	10 Yrs	15 Yrs
C-Stem AMT	Pinnacle	69	3025	1.3 (0.9, 1.8)	2.4 (1.8, 3.1)	3.1 (2.4, 4.1)	5.7 (3.6, 8.8)		

9. Berstock et al. A 5–8 Year Retrospective Follow-Up of the C-Stem AMT Femoral Component: Patient Reported Outcomes and Survivorship Analysis. *J. Arthroplasty* 2014 Sep; 29(9): 1753-7.
10. Flatoy et al. Triple taper stem design shows promising fixation and bone remodelling characteristics. Radiostereometric Analysis in a randomised controlled trial. *Bone Joint Journal* 2015 June; 97-B(6):755-61.
11. Stone M, et al. Wear Rate Analysis of Cemented Cross-Linked Polyethylene (Marathon) Acetabular Cups at Five Years. EHS 2014. Stockholm Poster No: P032.
12. Engh CA Jr, Stepniewski AS, Ginn SD, Beykirch SE, Sychterz-Terefenko CJ, Hopper RH Jr, Engh CA. "A randomized prospective evaluation of outcomes after total hip arthroplasty using cross-linked marathon and non-cross-linked Enduron polyethylene liners." *J Arthroplasty* 2006;21(6 Suppl 2):17-25.

This publication is not intended for distribution in the USA.

The third-party trademarks used herein are the trademarks of their respective owners.



Johnson & Johnson Medical Limited. Baird House, 4 Lower Gilmore Bank, Edinburgh, EH3 9QP, United Kingdom. Incorporated and registered in Scotland under company number SC132162.

DePuy Orthopaedics, Inc.
700 Orthopaedic Drive
Warsaw, IN 46582
USA
Tel: +1 (800) 366 8143
Fax: +1 (800) 669 2530

DePuy International Ltd
St Anthony's Road
Leeds LS11 8DT
England
Tel: +44 (0)113 270 0461

DePuy (Ireland)
Loughbeg
Ringaskiddy
Co. Cork
Ireland
Tel: +353 21 4914 000
Fax: +353 21 4914 199

DePuy Ireland UC
Loughbeg
Ringaskiddy
Co. Cork
Ireland
Tel: +353 21 4914 000
Fax: +353 21 4914 199



www.depuysynthes.com